


# Cardiovascular Sutures

- Our surgical sutures are made of high quality supplies from Europe.



## PACK

- ✓ Stability Assurance.
- ✓ Peel open.

## NEEDLE OF JAPANESE STEEL AISI 300 SERIES ORIGIN

- ✓ Silicone needle, which allows a single pass through the tissues.
- ✓ Designed to keep the edge after several passes.
- ✓ Needle of suitable size for each thread gauge.

## THREAD FROM GERMANY

- ✓ Uniform and homogeneous strand, which prevents tissue damage.
- ✓ High tensile strength.
- ✓ Flexible, facilitates handling and safety in performing the knots.
- ✓ Strand with appropriate length for each type of surgery.

## BENEFITS FOR THE PATIENT
















- Less tissue reaction and smooth passage through tissue due to the excellent edge and siliconized needles.
- Right concordance in the crimping that generates less trauma in the patient.

“ Suture made with advanced technology, which allows an excellent match between the thread and the needle. ”







# Polypropylene

Non absorbable monofilament synthetic suture.  
Color blue.

USP Size	Needle's shape and type	Needle size (*)	Length	
8/0	 3/8 circle micropoint 	6.5 mm	60 cm	
		8 mm		
7/0	 3/8 circle micropoint 	9.3 mm		
		10 mm		
6/0	 3/8 circle micropoint  	9.3 mm		75 cm
		10 mm		
5/0	 1/2 circle micropoint 	12 mm		
		15 mm		
4/0	 1/2 circle micropoint 	17 mm	75 cm - 90 cm	
		15 mm		
		20 mm		
		25 mm		
3/0	 1/2 circle micropoint 	30 mm		
		15 mm		
		20 mm		
2/0	 1/2 circle micropoint 	25 mm		
		20 mm		
		30 mm		





# Polyester / with and without PTFE Pledget

Non absorbable suture of synthetic, braided multifilament.  
Available in green - white / with and without PTFE Pledget.

USP Size	Needle's shape and type	Needle size (*)	Length
3/0	 1/2 circle micropoint 	15 mm	75 cm
		20 mm	
		25 mm	
2/0	 1/2 circle micropoint 	15 mm	75 cm
		20 mm	
		25 mm	

# Steel

Non absorbable suture of synthetic, monofilament.  
It has high permanent tensile strength.

USP Size	Needle's shape and type	Needle size (*)	Length
1	 1/2 circle round body cutting tip 	35 mm	4 x 45 cm
2		40 mm	
3		50 mm	
3	 1/2 circle round body cutting tip 	60 mm	
4			
5			

D/GDC-I-037 Rev. 01

**Recommendations:**

Open only sutures needed for the procedure. With gentle traction to straighten the sutures. Do not tighten, not rub, not abruptly pull the needles. Do not press or squeeze the threads of sutures with surgical instruments.

(\*) ± 2 mm.

**Warnings:**

Do not resterilize. Discard open unused sutures. Store at a temperature not exceeding of 27°C.