

# ENDO PLUS Sutures

It is the best suture to accelerate almost any ligation procedure.



The **ENDO PLUS** suture is ideal for ligation of pedicles used in laparoscopic procedures. It is composed of an 18" long ligature inserted in a plastic tube (narrow at one end and marked at the other).



The suture is formed in a ligation loop with a knot. Once the ligature is in place, Pull up to tighten the loop and secure the knot.

## PACK

- ✓ Stability Assurance.
- ✓ Peel open.

## NEEDLE OF STEEL AISI 302 SERIES


- ✓ Silicone needle, which allows a single pass through the tissues.
- ✓ Designed to keep the edge after several passes.
- ✓ Needle of suitable size for each thread gauge.

## THREAD

- ✓ Uniform and homogeneous strand, which prevents tissue damage.
- ✓ High tensile strength.
- ✓ Flexible, facilitates handling and safety in performing the knots.
- ✓ Strand with appropriate length for each type of surgery.

## BENEFITS FOR THE PATIENT

- Less tissue reaction and smooth passage through tissue due to the excellent edge and siliconized needles.
- Right concordance in the crimping that generates less trauma in the patient.



**FDA**  
Cleared

