

Microsurgery Sutures

Nervous Vascular

Reconstructive Procedures

ARTERIAL

- ✓ Term Suture - Lateral
- ✓ Term Suture - Terminal
- ✓ Lateral Arterial Suture - Terminal

VENOUS ANASTOMOSIS



NERVOUS MICROSURGERY

- ✓ Epineural
- ✓ Perineural
- ✓ Epi-Perineural

ARTERIAL INTERVENTION INJERSES

- ✓ Indo-Venous

PACK

- ✓ Stability Assurance.
- ✓ Peel open.

NEEDLE OF JAPANESE STEEL AISI 302 SERIES ORIGIN

- ✓ Silicone needle, which allows a single pass through the tissues.
- ✓ Designed to keep the edge after several passes.
- ✓ Needle of suitable size for each thread gauge.

THREAD FROM GERMANY

- ✓ Uniform and homogeneous strand, which prevents tissue damage.
- ✓ High tensile strength.
- ✓ Flexible, facilitates handling and safety in performing the knots.
- ✓ Strand with appropriate length for each type of surgery.



BENEFITS FOR THE PATIENT







- Less tissue reaction and smooth passage through tissue due to the excellent edge and siliconized needles.
- Right concordance in the crimping that generates less trauma in the patient.



Black Nylon

Non absorbable monofilament suture structure

USP Size	Needle's shape and type	Needle size (*)	Length
9/0, 10/0	 3/8 circle micropoint 	4.0 mm, 5.0 mm, 5.5 mm	30 cm

USP Size	Needle's shape and type	Needle size (*)	Length
10/0	 3/8 circle micropoint 	4.0 mm,	15 cm
10/0 9/0	 3/8 circle micropoint 	5.0 mm, 5.5 mm	15 cm
8/0	 3/8 circle micropoint 	6.5 mm	15 cm

“ Suture made with advanced technology, which allows an excellent match between the thread and the needle. ”



D/GDC-043; Rev. 01

Recommendations:

Open only sutures needed for the procedure. With gentle traction to straighten the sutures. Do not tighten, not rub, not abruptly pull the needles. Do not press or squeeze the threads of sutures with surgical instruments.
(*) ± 2 mm.

Warnings:

Do not resterilize. Discard open unused sutures. Store at a temperature not exceeding of 27 °C.