

# Ophthalmic Sutures

- Manufactured with high quality supplies.
- Specially designed for eye diseases and others surgical procedures.



## PACK

- ✓ Stability Assurance.
- ✓ Peel open.

## NEEDLE OF JAPANESE STEEL AISI 300 SERIES ORIGIN

- ✓ Silicone needle, which allows a single pass through the tissues.
- ✓ Designed to keep the edge after several passes.
- ✓ Needle of suitable size for each thread gauge.

## THREAD FROM GERMANY

- ✓ Uniform and homogeneous strand, which prevents tissue damage.
- ✓ High tensile strength.
- ✓ Flexible, facilitates handling and safety in performing the knots.
- ✓ Strand with appropriate length for each type of surgery.

## BENEFITS FOR THE PATIENT





- Less tissue reaction and smooth passage through tissue due to the excellent edge and siliconized needles.
- Right concordance in the crimping that generates less trauma in the patient.

“ Suture made with advanced technology, which allows an excellent match between the thread and the needle. ”





## Polyglycolic Acid

Multifilament absorbable suture structure

USP Size	Needle's shape and type	Needle size (*)	Length
6/0, 7/0, 8/0	 3/8 spatulate circle 	6.2 mm	30 / 45 cm
6/0	 1/4 spatulate circle 	8.0 mm	30 / 45 cm





## Black Nylon

Non absorbable monofilament suture structure

USP Size	Needle's shape and type	Needle size (*)	Length
9/0, 10/0	 3/8 spatulate circle 	6.0 mm, 6.2 mm, 6.4 mm	30 / 40 cm
8/0			30 cm

## Black Silk

Non absorbable suture braided structure

USP Size	Needle's shape and type	Needle size (*)	Length
6/0	 1/4 spatulate circle 	6.2 mm, 8.0 mm	45 cm
7/0, 8/0	 3/8 spatulate circle 	6.2 mm	45 cm

## Virgin Silk

Non absorbable suture twisted structure

USP Size	Needle's shape and type	Needle size (*)	Length
7/0, 8/0	 3/8 spatulate circle 	6.0 mm, 6.2 mm	45 cm

### FREQUENTLY USES:

In corneal surgery, cataract, glaucoma, retina, vitrectomy, lio fixation, strabismus, plastic surgery.

Certifications which feature our products:



#### Recommendations:

Open only sutures needed for the procedure. With gentle traction to straighten the sutures. Do not tighten, not rub, not abruptly pull the needles. Do not press or squeeze the threads of sutures with surgical instruments.  
(\*) ± 2 mm.

#### Warnings:

Do not resterilize. Discard open unused sutures. Store at a temperature not exceeding of 27°C.