

# SURGICAL SUTURES

**SUTUMED** 

Medical Devices

U.S. Surgical Suture Manufacturer

## We manufacture surgical sutures with highest quality

Our product has been manufactured according to the controls established by a quality management system that meets the requirements of international certifications.

### Our surgical sutures fulfill most important benefits for the surgeon.

- ★ Sterile.
- ★ Easy to handle.
- ★ Without propensity for bacterial growth.
- ★ Able to resist when knotted without fraying or cutting.
- ★ Minimal tissue reaction after finished its purpose.


**FDA** **CE** **GMP**  
CLEARED

**ISO** **ISO**  
9001 13485




### Needle's Shape

Thinner, stronger and better concordance between the needle and the thread




### Thread Benefits

Concentric crimping that produce a edge minor to pass through tissue. This allows the surgeon greater control of surgery and less tissue trauma to the patient.



### Primary Label

- 1- Metric
- 2- Strand length
- 3- Shape of the needle
- 4- Sterilization
- 5- Expiration date
- 6- Lot
- 7- Absorbable / non absorbable
- 8- Material
- 9- Use



- 10- Needle's shape and size
- 11- USP size




**SUTUMED**   
Medical Devices


SUTUMED CORP


Principal Executive office

 Florida, US

 9280 College Pkwy, Unit# 5-6 Fort Myers, FL 33919

 (239) 590 8484

 sutumed.us

 info@sutumed.us

## Absorbable Sutures

FDA CE GMP ISO ISO  
CLEARED 2121 9001 13485

### Suture Characteristics


**SutuCRYL / POLYGLACTIN**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**SutuCRYL POLYGLACTIN**  
Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: By hydrolysis  
COLOR: Violet, undyed  
USP SIZE: 8-0 to 2.

STRUCTURE: Braided Multifilament




**SutuCRYL MONO / POLYGLYCAPRONE**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**SutuCRYL MONO POLYGLYCAPRONE**  
Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: By hydrolysis  
COLOR: Violet, undyed  
USP SIZE: 5-0 to 1.

STRUCTURE: Monofilament




**CHROMIC GUT**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**CHROMIC GUT CATGUT CRÓMICO**  
Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: By phagocytosis  
COLOR: Bright brown, dark brown  
USP SIZE: 6-0 to 3.

STRUCTURE: Twisted Multifilament



### Suture Characteristics


**PDO II / POLYDIOXANONE**

2-0 (3 Metric) 28"(70 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**PDO II POLYDIOXANONE POLIDIOXANONA**  
Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: By hydrolysis  
COLOR: Violet  
USP SIZE: 7-0 to 2.

STRUCTURE: Monofilament




**SutuSORB / POLYGLYCOLIC ACID**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**SutuSORB POLYGLYCOLIC ACID ÁCIDO POLIGLICÓLICO**  
Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: By hydrolysis  
COLOR: Violet  
USP SIZE: 8-0 to 2.

STRUCTURE: Braided Multifilament




**PLAIN GUT**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**PLAIN GUT CATGUT SIMPLE**  
Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: By phagocytosis  
COLOR: Cream, bright yellow  
USP SIZE: 6-0 to 3.

STRUCTURE: Twisted Multifilament



## Non-absorbable Sutures

### Suture Characteristics


**SutuLON / MONOFILAMENT NYLON**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**SutuLON MONOFILAMENT NYLON NYLON MONOFILAMENTO**  
No Abs. / Surgical Suture  
Sutura Quirúrgica USP  
BLACK / NEGRO

TYPE OF ABSORPTION: Non-absorbable  
COLOR: Black  
USP SIZE: 10-0 to 2.

STRUCTURE: Monofilament




**SutuSILK / VIRGIN SILK**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**SutuSILK VIRGIN SILK SEDA VIRGEN**  
No Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: Non-absorbable  
COLOR: Blue, black  
USP SIZE: 8-0 to 7-0.

STRUCTURE: Twisted Multifilament




**SutuSTEEL / STEEL**

5 (0.5 Metric) 18"(45 cm) ST8755S  
STBV1754  
Taper

**SutuSTEEL STEEL ACERO**  
No Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: Non-absorbable  
COLOR: Silver  
USP SIZE: 1 to 5.

STRUCTURE: Monofilament and Multifilament



### Suture Characteristics


**SutuSILK / BLACK BRAIDED SILK**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**SutuSILK BLACK BRAIDED SILK SEDA NEGRA TRENZADA**  
No Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: Non-absorbable  
COLOR: Black, White  
USP SIZE: 7-0 to 2.

STRUCTURE: Braided Multifilament




**SutuLENE / POLYPROPYLENE**

7-0 (0.5 Metric) 30"(75 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**SutuLENE POLYPROPYLENE POLIPROPILENO**  
No Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: Non-absorbable  
COLOR: Blue  
USP SIZE: 10-0 to 1.

STRUCTURE: Monofilament



**SutuESTER / POLYESTER**

2-0 (3 Metric) 18"(45 cm) ST8755S  
STBV1754  
15 mm 1/2 Taper Point

**SutuESTER POLYESTER POLIESTER**  
No Abs. / Surgical Suture  
Sutura Quirúrgica USP

TYPE OF ABSORPTION: Non-absorbable  
COLOR: Green, White  
USP SIZE: 6-0 a 5.

STRUCTURE: Braided Multifilament

